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	Standard soda lime	New generation absorbent containing traces of alkali-hydroxide	New generation Containing no alkali-hydroxide
Alkali Hydroxide (normally Sodium Hydroxide)	3 %	0.5 % to 3 %	Nil
Humectant (to retain moisture)	Nil	0.5 % to 8 %	0.5 % to 8 %
pH sensitive indicator	Trace	Trace	Trace
Calcium Hydroxide	>96 %	>90 %	> 90 %

Note, these figures represent the dry weight constituents. Water makes up 13 % to 17 % of the finished product.

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09/16/1996: While undergoing surgery, pt developed high COHb level of 58%. The anesthesia machine had not been used for a period of 36 to 48 hours~~~Ongoing investigation has resulted in the conclusion that the event was a result of the reaction between dry baralyme and desflurane.

10/30/1996: While undergoing surgery, patient developed high COHb level of 49%. The machine had not been used for about 48 hrs. ~~~.It is considered that the high CO level was the result of an interaction between baralyme and desflurane.

09/25/2003: The absorber of the anesthesia machine caught on fire. The machine was not in use at the time, however the oxygen and sevoflurane had been left on~~~ There were two baralyme canisters in the anesthesia machine.

Ref 3: Reports found on the FDA Medical Device Database

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JA Baum, JH Woehlck
chemical reaction Interaction of
Inhalation anaesthetics with CO₂ absorbents. Best
Practice & research Clinical Anaesthesiology. Vol.
17. No. 1. Pp. 63-76, 2003


Figure 2. Utilization time of different carbon dioxide absorbents, determined in clinical practice during routine performance of minimal and closed-system anaesthesia. Legend downwards: *n* = number of the respective measurement runs; mean utilization time (hours); range of utilization times; MFT = percentage of time using extremely low fresh gas flow rates between 0.5 and 0.25 l/min. The wide range of utilization times of Amsorb reflects the differences in absorption capacity of three different batches of this absorbent.

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